

USER INSTRUCTIONS

CLEANOR 304 ELECTROPOLISHING FOR GOLD					
Colour : —	Use : Bath	Colour coordinates : L: —	a: —	b: —	

Product description

CLEANOR 304 is a special electropolishing process ready to use for gold. It is completely free of cyanides and heavy metals and can polish effectively, ensuring minimum treatment loss compared to other similar products on the market.

Operating conditions

Parameter	Average suggested value	
Treatment time	500 seconds	
Temperature	60 °C	
Voltage	8 V	
Current density	_	
Cathodic efficiency	_	
Deposition speed	_	
Anodes	Stainless Steel 316 or graphite	
Agitation	Necessary	
Ventilation	Necessary	

Maintenance of plating bath

The consumption of solution is directly proportional to the amount of polished goods. We recommend changing it when the solution is dirty or when there are residues at the bottom.

Safety Information

Please refer to the safety data sheets for any information regarding safety and disposal (MSDS).



USER INSTRUCTIONS

CLEANOR 304 ELECTROPOLISHING FOR GOLD				
Colour : —	Use : —	Colour coordinates : L: —	a: —	b: —

NOTE

The electropolishing process works with the opposite polarity to that of a traditional galvanic bath. In the electropolishing process, the objects to be treated are connected to the POSITIVE pole of the voltage rectifier, while the anodes are connected to the NEGATIVE pole.

Notes

Instructions indicated in Our Process technical data sheets are the result of attentive checks and have been written as a guideline. They represent the best of our knowledge and refer to a normal use of our products. Such products are guaranteed in their quality up to the delivery: we cannot guarantee for the correct use of them, as this action in not under our direct control (addiction or correction to plating baths are made by customer's personnel). As a consequence, process instructions have not to be intended as a guarantee of the final result in the use of our products.

Related products					
Cod.Articles	Description	Colour	Packaging		
EB010201	CLEANOR 304 electropolishing for gold (ready to use)	_	25 L		